



# Image j/Fiji for beginners


Mathieu Fallet

June 2017



# Plan


- Introduction
- Help menu
- Task bar
- File menu
- Image menu
- Process menu
- Analyse menu
- Plugins menu




# Introduction

- Fiji is Image J with plugins inside
- It is actually based on java 8
- Fiji integrate the new Image J2 which have better functionalities (plugins update, script editor, better reading and writing image format)
- Plugins or macro should be downloaded in the plugin directory of Fiji
- Lot of information here :  
<https://imagej.nih.gov/ij/docs/guide/146-29.html>


# Help menu : to begin



- Documentation : [Help/Documentation](#)
- Plugins: [Help/Plugins](#)
- Update : [Help/update Image J and update Fiji](#)
- Use Ctrl 1 to search function
- Options : [Edit/Options/Memory and Threads and Image J2](#)




# Taskbar




To see the option double clic or clic right



- Shape : to define region and crop
- Line : to draw line
- Points : to count
- Magic Wand : to define object
- Annotation : to draw annotation (ctrl D to draw it)
- Loop and hand
- Color Picker : to define foreground and background color
- Menu




# File menu



- Drag and drop your images in the taskbar
- Import multiple TIFF : [File/Import/Image sequence](#)
- Plugin Bio-formats
- Save As (avi...) no compression use VirtualDub to compress in mpeg4
- Make a new image : [File/New Image](#)
- Revert : [File/Revert](#)
- Open : [File/Open Samples](#), [Open recent](#)




# Edit menu





- Fill a space : [Edit/Fill \(ctrl D\)](#)
- Draw a text : [Edit/Draw \(ctrl F\)](#)
- Clear : [Edit/Clear](#) and [Clear outside](#)
- Make selection : [Edit/Selection/Select all](#) or [none](#) or [restore selection](#) or [Create Selection](#) (on binary mask) or [Specify](#) (to crop on specified number of pixel)
- Options : [Edit/Options/Memory and Threads](#) and [Image J2](#)
- Invert : [Edit/invert](#)

# Image menu




- Check calibration : [Image/show Info and properties](#)
- Handle colors channels : [Image/colors/channel tool](#),  
[image/color/split channels](#)
- Change LUT : [Image/lookup Tables](#)
- Change contrast : [Image/adjust/Brightness&contrast](#) or  
[image/adjust/color balance](#)
- Change type images : [Image/type](#)
- Handle stack images : [Image/stacks/tools](#)
- Crop : [Image/Crop](#)









# Process menu



- Median filters : [Process/Filters/median](#)
- Morphology filters : [Process/binary \(options\)](#)
- FFT : [Process/FFT \(Bandpass Filter, FD math\)](#)
- Find spot : [Process/Find Maxima](#)




Points at maxima (Multi-point selection)






Segmented Particles (ROIs obtained with  
[Analyze>Analyze Particles...](#))


# Process menu






- Mathematical operation : [Process/Math](#)
- Operation between images : [Process/Image calculator](#)
- Subtract background : [Process/Subtract background](#)




# Analyse menu




- Measure : Analyse/measure and set measurements
- Analyse : Analyse/Analyse particles...
- Histogram : Analyse/Histogram
- Profile : Analyse/Plot Profile
- Tools : Analyse/tools/Scale bar and ROI manager



# Plugin and macro



- Macro : [Plugins/Macros/Record...](#)
- Shortcut : [Plugin/Shortcut](#)
- Control panel : [Plugin/Utilities/Control panel...](#)
- List of plugins: try Jacop for colocalisation
- Possibility to add icon on the taskbar



# Shortcut


## Interesting shortcuts: ctrl

a	Select All
A	Select None
b	Add Selection...
C	Brightness/Contrast...
c	Copy
d	Draw
D	Duplicate...
E	Restore Selection
e	Scale...
f	Fill
F	Flatten
G	Capture Screen
h	Histogram
H	Orthogonal Views
I	Invert
i	Show Info...
K	Color Picker...

k	Plot Profile
<b>I</b>	<b>Find Commands...</b>
M	Install...
m	Measure
N	Text Window
o	Open...
P	Properties...
r	Revert
R	Repeat last command
s	Save
S	Smooth
t	Add to Manager
T	Threshold...
U	Control Panel...
v	Paste
V	System Clipboard
w	Close
X	Crop
x	Cut
Z	Channels Tool...
z	Undo


Open the image “Thymus” :

- Change the contrast channels
- Change the LUT of the blue channels in cyan
- Make the Z projection, make a median filter
- Draw a scale, check the calibration



Open the image Noyaux :

- Split the channel
- Make the threshold on the blue channel and create the binary mask
- Do the segmentation and count the objects
- Do a macro recording the different steps



# Practice

Open the image DICand correct the background by creating a shading image.

Open the image EV2NK and apply a median filter. Calculate the Pearson coefficient with and without median filter.

Open the image ABCA1 and quantify the membrane staining. How to do it ?

